

Digital Rocker

Digital display, large platform, speed control, and rapid angle adjustment for reliable sample agitation



The Thermo Scientific™ Digital Rocker is designed for ELISA, protein synthesis, gel staining and de-staining, antibody staining, and other time-sensitive applications.

Key Features

- Digital speed and time displays
- Adjustable rocking angle: 1° to 15°
- Speed range: 1 to 70rpm
- Speed control:
 - ≤50rpm: ±1rpm
 - >50rpm: ±2rpm
- Timer range: continuous or 1 min. to 99 hr. 59 min.
- Capacity: 6 x 96-well plates
- Maximum load: 4.5kg at ≤50rpm




Includes Dimpled Rubber Mat (88882102) and Rubber Strips (88882103, 6/Pack). Optional non-slip rubber mat tray is designed for carrying plates, dishes, culture bottles, and bags. Warranty: 2 years | Certifications: cCSAus, CE

Digital Rocker shown with flasks on top of a Thermo Scientific Dimpled Rubber Mat. Flasks sold separately.

| Product | Electrical Requirements | Tray Dimensions D x W x H | Overall Dimensions D x W x H | Weight | Plug Type | Catalog Number |
|----------------|-------------------------|--|---|----------------------|----------------|-----------------|
| Digital Rocker | 100 to 240V, 50/60Hz | 32.2 x 28.2 x 1.9cm (12.7 x 11.1 x 0.7 in.) | 39 x 29.2 x 20.6cm (15.4 x 11.5 x 8.1 in.) | 13.5kg (29.8 lb.) | 120V USUP | 88882001 |
| Digital Rocker | | | | | 230V EU/UK/CHN | 88882002 |

Accessories

| Product and Description | Quantity | Catalog Number |
|--|----------|-----------------|
| Rubber Mat, 32.2 x 28.2cm (12.7 x 11.1 in) | Each | 88882101 |
| Dimpled Rubber Mat, 32.2 x 28.2cm (12.7 x 11.1 in) | Each | 88882102 |
| Rubber Strips | 6/pack | 88882103 |
| Power Adapter with US plug | Each | 88870126 |
| Power Adapter with AUS/CN plug | Each | 88870127 |
| Power Adapter with EU plug | Each | 88870128 |
| Power Adapter with UK plug | Each | 88870129 |



Thermo Scientific Rubber Strips (top), Rubber Mat (bottom left), and Dimpled Rubber Mat (bottom right).